

SPUTUM COLLECTION FOR BK MICROBIOLOGICAL ANALYSIS

Procedure no.: MIC-IPA-011

Version: 003

Status: Active

1. Sample collection

- Do not rinse your mouth. Spit into the sterile wide-neck container provided by the laboratory.
- Make sure not to spit saliva or nasal secretions. This will make the collection unusable.
- Close the container tightly to prevent leaking, making sure not to touch the inside of the container.
- On the container, write your name, file number and the date and time of the collection. Then place the container in a Ziploc bag.
- Ideally, do 3 collections on the same day at intervals of at least one hour.
- The best time to collect the samples is in the morning before breakfast or within 24 hours following a bronchoscopy.

2. Storage

- This type of sample **is very difficult to store.**
- Keep the sample **in the refrigerator** until transporting it to the laboratory reception (SBA-31), **as soon as possible.**

3. Transportation

- Bring the sample, accompanied by the requisition form duly completed, to the laboratory reception (SBA-31), **as soon as possible** (ideally within two hours of collection).
- Maximum timeframe: 12 hours after the collection, if refrigerated.
- For better storage, keep the sample **cool** during transportation, in a cooler with an ice pack.

4. For more information

- Contact the laboratory reception at 450-759-8222, ext. 2606 (business days, 8 a.m. to 4 p.m.).

Martine Lavergne, MD
Microbiologist-infectious disease
specialist

Jean-François Boisvert, MD
Microbiologist-infectious disease
specialist

Soraya Boukhoudmi, MD
Microbiologist-infectious disease
specialist

Anaïs Lauzon-Laurin, MD
Microbiologist-infectious disease
specialist