

STOOL CULTURE COLLECTION

Procedure no.: MIC-IPA-004

Version: 007

Status: Active

1. Sample collection

- It is important to not let the stool touch the water in the toilet bowl and to not use toilet paper to collect it.
- Use a basin or a clean large container to collect the stool.
- Collect 3 to 6 g of solid stool (about the size of a grape) or 5 mL of liquid stool and place it in the Cary-Blair transport medium provided by the laboratory (**do not pass the line on the container**).
- Close the container tightly to make sure that it does not leak.
- Shake the container vigorously for 10 seconds to make sure that the stool is broken up into small pieces and well combined with the transport medium.
- Check the appropriate box on the container to indicate if the stool is solid, soft, liquid or contains blood.
- On the container, write your name, file number and the date and time of the collection. Then place the container in a Ziploc bag.
- Only one collection per day is accepted. Other collections taken on the same day will be rejected.

2. Storage

- Keep the sample at **room temperature (20-25°C)** until transportation. Do not refrigerate.

3. Transportation

- Bring the sample, accompanied by the requisition form duly completed, to the laboratory reception (SBA-31), **as soon as possible**.
- Maximum timeframe: 24 hours before collection.

4. For more information

- Contact the laboratory reception at 450-759-8222, ext. 2606 (business days, 8 a.m. to 4 p.m.).

Martine Lavergne, MD
Microbiologist-infectious disease
specialist

Jean-François Boisvert, MD
Microbiologist-infectious disease
specialist

Soraya Boukhoudmi, MD
Microbiologist-infectious disease
specialist

Anaïs Lauzon-Laurin, MD
Microbiologist-infectious disease specialist

